

Simply a good feeling. System Hygiene from Dürr Dental.

4 colours, 4 areas - simply logical, simply safe



It is best when you can rely on a system

In 1965 Dürr Dental launched Orotol[®], the first disinfectant for suction systems. This marked the beginning of the continuous development and refinement of our expertise in the hygiene area. Today, dentists throughout the world trust Dürr Dental System Hygiene to ensure comprehensive protection for their surgeries.


Consulting

- Training courses, staff training sessions and on-site consulting
- Expert knowledge from qualified trainers
- Solution-oriented practice tips tailored to individual dental surgeries
- Support for development of a quality management system



Service

- Professional and solution-oriented product advice
- Up-to-date online information



System Hygiene from Dürr Dental




Quality

- High production standards - Made in Germany
- Defined and controlled work processes ensure safety
- Decades of experience in the manufacture of disinfectants and cleaning agents



In addition to Dürr Dental System Hygiene products, we also offer a comprehensive service package to support your practice team in their daily work. We stand by your team with our on-site consulting service. Our service department will help you with questions about our products. And outstanding quality as the product of decades of experience goes hand-in-hand with frequent innovations to meet the needs of dental practices.



User orientation

- Guidance thanks to the System Hygiene colour system
- Distinctive and easy-to-store 2.5 l bottle
- Ergonomic bottle shape makes for easy handling
- Application aids for perfect dosing



Hygiene and technology

- Over 80 years of experience in medical technology
- Unbeatable: technical solutions in combination with the cleaning and disinfecting solutions from Dürr Dental
- Gentle and effective products shine with easy handling and speed of use



4 colours - 4 areas

- Blue for instruments **Page 08**
- Green for surfaces **Page 12**
- Pink for skin and hands **Page 16**
- Yellow for special areas **Page 18**

It is best when a single system offers a complete solution

When it comes to infection protection, disinfection and cleaning, nothing should be left to chance. This is why we have carried out countless tests to demonstrate that all our preparations meet the key demands of the day-to-day surgery environment.

Safe

- Comprehensive spectrum of activity: bactericidal, tuberculocidal, fungicidal, virucidal incl. coronaviruses
- Certified, reliable effect
- VAH-listed and EN-compliant disinfectants, without exception
- High disinfecting and cleaning power

Our recommendation for your dental surgery:

HD 410
Hand disinfection

HD 435
Washing lotion

FD 366 sensitive
Disinfection of sensitive surfaces

FD 360
Cleaning and care of vinyl upholstery

FD 366 sensitive wipes
or sensitive top wipes

FD 322 or FD 333
Quick-acting
disinfection

FD 333
forte premium wipes

FD 333 top wipes

Gentle


- Superior combination of active agents
- Particularly gentle and kind to skin
- Suitable for use on most surfaces and recommended by leading equipment manufacturers
- Clinically tested
- Bio-degradable

Convenient

- Easy and efficient to use
- Convenient, standard container sizes
- Cleaning according to plan
- Dosing caps for practical mixing
- Highly cost-effective to use

Fast

- Short reaction times
- Quick, residue-free drying
- Practical application



MD 550
Spittoon bowl cleaner



Orotol® plus
Disinfection of Suction Systems

MD 555 cleaner
Special detergent for
suction systems



ID 212 or ID 213
Instrument disinfection



ID 220
Bur disinfection



Dürr Dental Hygiene technology – for maximum efficiency in the hygiene cycle

For more than 80 years, Dürr Dental has been synonymous with trend-setting developments in medical technology – Made in Germany. With this tradition and the goal of systematically perfecting our range of hygiene technology products, today we can offer a variety of innovative technical solutions


Use

Manual cleaning and disinfection

Machine reprocessing in the washer-disinfector

Function test/check and instrument care

Sterile packaging with Hygopac View and Hygofol



so that day-to-day hygiene and preparatory measures in dental surgeries can be performed most efficiently and with the best results. In Dürr Dental System Hygiene, the individual components work together to create an unbeatable team.

Steam sterilisation in
the Hygoclave 40/50

Documentation

Storage



Instruments



Disinfection and cleaning Instruments and drills

■ ID 212 Instrument disinfection

Concentrate for cleaning and disinfecting instruments with excellent material compatibility



- Bactericidal, tuberculocidal in an ultrasonic device, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), as well as non-enveloped adenoviruses and noroviruses
- Suitable for general, surgical and rotary instruments
- Very short reaction time
- Outstanding material compatibility with instruments containing rubber and silicon

Active ingredients: Quaternary ammonium compounds

Available as: 2.5 l

Our recommendation: Hygobox

■ ID 212 forte Instrument disinfection

Concentrate for cleaning and disinfecting general, surgical and rotating instruments



- Bactericidal, tuberculocidal in an ultrasonic device, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), as well as non-enveloped adenoviruses and noroviruses
- Suitable for general, surgical and rotary instruments
- High corrosion protection and long working time of 7 days
- Fast disinfection and cleaning

Active ingredients: Alkylamines, quaternary ammonium compounds

Available as: 2.5 l

Our recommendation: Hygobox

■ ID 213 Instrument disinfection

Concentrate for cleaning and disinfecting general, surgical and rotating instruments with a broad range of effectiveness



- Bactericidal, tuberculocidal, mycobactericidal in an ultrasonic device, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV), noroviruses and polyomaviruses SV40
- Very economical thanks to a long working time of 14 days
- Excellent material compatibility thanks to corrosion inhibitors

Active ingredients: Alkylamines, quaternary ammonium compounds

Available as: 2.5 l

Our recommendation: Hygobox

■ ID 215 Enzymatic instrument cleaner

Concentrate for manual, enzymatic instrument cleaning prior to disinfection



- Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV)
- With proven personnel protection
- Exceptional cleaning efficiency - dissolves even the most stubborn contamination
- Good material compatibility

Active ingredients: Quaternary ammonium compounds, guanidine compounds

Available as: 2.5 l

Our recommendation: Hygobox

■ ID 220 Bur disinfection

Fully virucidal, ready-to-use solution for cleaning and disinfecting rotating instruments



Tested in accordance with the DVV/RKI Directive as well as EN 14476 and EN 17111

- Bactericidal, mycobactericidal, tuberculocidal, yeasticidal, fungicidal, virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses, as well as non-enveloped viruses such as adenoviruses, polyomaviruses SV40, noroviruses and polioviruses)
- Exceptional cleaning power, especially designed for the microstructures in drills and cutters
- Outstanding corrosion protection and long working time of 7 days
- Very short reaction times

Active ingredients: Alcohols in alkaline medium

Available as: 2.5 l



Devices and accessories

■ Hygobox

Sturdy and robust transport and disinfection box with a 3 litre capacity



- Ergonomic handles for safe transport
- Tilt-resistant drip-off position for the integrated cartridge
- Specially shaped slits prevent the instruments from becoming entangled

Dimensions (HxWxD): 13 x 32.5 x 21 cm

Colours: Anthracite-coloured lid with blue or yellow cartridge, lid available in 5 different colours and with a white cartridge

Our recommendation: ID 212, ID 212 forte, ID 213, ID 215, MD 520, MD 535

■ Hygopac

Rotary sealing device for transparent sterile packaging



- Optimal pressing pressure for consistent sealing seam quality
- Assembly-line principle for fast and economical operation
- Continuous temperature regulation for optimal seal seam

Dimensions (HxWxD): 16.0 x 43.5 x 13.5 cm

Our recommendation: Hygofol bags, Hygofol rolls, instrument table, Hygofol station

■ Hygopac View

Validatable rotary sealing device for standard-compliant sealing of sterilised packaging in accordance with DIN EN ISO 11607-2.



- Exceptionally compact design with minimal space requirements
- Intuitive user interface via touch display
- Process can be documented by saving to an SD card and via the network connection

Dimensions (HxWxD): 17.0 x 37.0 x 14.0 cm

Our recommendation: Hygofol bags, Hygofol rolls, instrument table, Hygofol station



Devices and accessories

Hygoclave 40/50

With the Hygoclave 40 and 50 series, Dürr Dental offers dental practices a professional Class B solution that combines impressive performance with a well-thought-out operating concept – for maximum efficiency and exceptional user friendliness designed to cope effortlessly with tough day-to-day working environments.



- High-performance sterilisation and drying system
- Excellent usability thanks to:
 - Intuitive operating concept
 - High-resolution 5" touch screen
 - Automatic door locking mechanism
- Flexible fresh water supply and integrated fresh water quality control
- Low water and energy consumption
- Everything in view and documented thanks to connection to the practice network and dental documentation programmes*

Dimensions:

Hygoclave 40 (HxWxD): 49.3 x 48.4 x 65.4 cm

Hygoclave 40 Plus (HxWxD): 49.3 x 48.4 x 65.4 cm

Hygoclave 50: 49.3 x 48.4 x 64.2 mm

Hygoclave 50 Plus: 49.3 x 48.4 x 64.2 mm

*Only available for the Hygoclave 50 series

Hygotest

The Helix test phantom with 250 test strips for a regular, quick and simple Helix test.



- In accordance with EN 867-5 and ISO 11140-1
- PCD (Process Challenge Device) and 250 indicators
- Can be used as a batch control system and for a Bowie-Dick test
- Starter set with 50 indicators included with the initial autoclave delivery

Hygofol

Hygofol sterile product packaging offers maximum reliability with its multi-layer laminated foil and highly tear-resistant sterilisation paper



- High puncture resistance and excellent mechanical strength
- High sealing seam tensile strength
- Hygofol as practical, prefabricated bag product in 3 sizes
- Hygofol for individual packaging as roll product in 5 sizes

Dimensions:

Hygofol bag (WxL): 7.5 cm, 10 cm or 15 cm x 30 cm each

Hygofol rolls (WxL): 5 cm, 7.5 cm, 10 cm, 15 cm or 25 cm x 100 m each



■ Hygoprint Easy

The optionally available printer offers the possibility to print cycle protocols as well as barcode labels. These contain all the required information on sterilisation cycles completed.



- Label size in accordance with guidelines - can be applied outside the sealing seam
- Bar code for scanning in digital data management
- Label specifications: device (designation and serial number), program, batch, date, user (approval by), expiration date

Dimensions :

Hygoprint Easy (HxWxD): 16.9 x 12.5 x 22.1 cm

Hygoprint Easy labels: (HxW): 2.5 cm x 5.1 cm

Our recommendation: Hygoprint Easy labels and paper rolls

■ Hygoprint

System for labelling sterile product packaging.



- Quick and simple documentation of the sterilisation process
- Complete labelling with all relevant information
- Very lightweight and handy

Our recommendation: Single or double-sided adhesive labels, documentation block with 100 sheets, replacement ink roller

■ Hygodem 90

Hygodem 90 delivers high quality in a compact unit. The ion exchanger unit filters tap water in three cartridges and can be directly connected to the Hygoclave 50, requiring little space.



- Ion exchange unit
- Convenient, thanks to 3 ion exchanger cartridges
- Simple maintenance, thanks to relief valve
- Simple installation, thanks to quick connector
- Low space requirements
- Affordable purchase price

Our recommendation: Cartridge set for Hygodem 90

■ Hygopure 90

When operated at maximum performance in the harsh day-to-day practice environment, Hygoclave needs a powerful and reliable partner. With its large capacity, Hygopure 90 enables removal of large quantities (several litres) of DI water at once.



- Reverse osmosis unit incl. removal tap and storage tank
- Multiple filter stages + high-performance membrane with a service valve for checking effectiveness
- Quick connector makes for simple installation
- High capacity
- Multiple litres of DI water can be removed at once
- The eco module permits water savings of up to 70%

Our recommendation: Cartridge set for Hygopure 90, DF50 membrane for Hygopure 90, stainless steel removal tap for Hygopure 90, HMEC1 key valve measurement device, LD1 Aqua-Stop valve

Surfaces



Disinfection and cleaning for surface disinfection

■ FD 300 Surface disinfection

Fully virucidal and alcohol-free concentrate for cleaning and disinfection of surfaces of non-invasive medical devices and fixtures



- Bactericidal, yeasticidal, tuberculocidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), virucidal (against enveloped viruses and non-enveloped viruses).
- Especially fast reaction time
- Economical application: Up to 250 l of ready-to-use solution at 1 %
- Good material compatibility

Active ingredients: Alkylamines, quaternary ammonium compounds

Available as: 2.5 l

Our recommendation: FD multi wipes, FD multi wipes compact

■ FD 312 Surface disinfection *

Alcohol-free concentrate for cleaning and disinfection of surfaces of non-invasive medical devices, fixtures and floors



- Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses).
- Economical application: Up to 250 l of ready-to-use solution at 1 %
- Good material compatibility
- Pleasant scent

Active ingredients: Quaternary ammonium compounds

Available as: 2.5 l

Our recommendation: FD multi wipes, FD multi wipes compact

Quick-acting disinfection

■ FD 322 Quick-acting disinfection

Fast-acting, ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40.
- Very short reaction time
- Dries without residue
- Good material compatibility

Active ingredients: Alcohols

FD 322 available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l

FD 322 perfume-free available as: 2.5 l, 10 l

Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact



■ FD 333 Quick-acting disinfection

Fast-acting, ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses).
- Good material compatibility
- Fresh lemon scent
- Quick drying

Active ingredients: Alcohols, quaternary ammonium compounds

Available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l

Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

■ FD 333 forte Quick-acting disinfection

Fully virucidal ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, mycobactericidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), virucidal (enveloped and non-enveloped viruses).
- Good material compatibility
- Fresh lemon scent
- Quick-acting disinfection effect

Active ingredients: Alcohols, quaternary ammonium compounds

Available as: 1 l, 2.5 l, 10 l

Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

■ FD 366 sensitive Disinfection of sensitive surfaces

Alcohol-reduced ready-to-use solution for the quick disinfection of sensitive surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses.
- Quick-acting disinfection effect
- Ideal for treatment units in combination with FD 360 Cleaning and care of vinyl upholstery
- Certified material compatibility – recommended by leading equipment manufacturers

Active ingredients: Alcohols, quaternary ammonium compounds

Available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l

Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

Care

■ FD 360 Cleaning and care of vinyl upholstery

Cleaning fluid for deep cleaning and care of vinyl upholstery



- Cleaning agent, gentle on materials, combined with a care product for long-term value-retention
- Intensive cleaning effect - removes stains of all kinds, rime and discolouration
- Certified material compatibility – recommended by leading equipment manufacturers
- Ideal for treatment units in combination with FD 366 sensitive

Ingredients: Special surfactants, silicon compounds and care formula based on avocado oil

Available as: Set with 500 ml spray bottle, 1 FD polishing wipe and 9 special sponges

Our recommendation: FD 366 sensitive, FD 366 sensitive wipes, FD 366 sensitive top wipes, FD polishing wipe

Surface disinfection wipes

Alcohol-free disinfection wipes

■ FD 312 wet wipes* Surface disinfection

Alcohol-free disinfection and cleaning wipes for floors, surgery inventory etc.



- Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses)
- High yielding surface coverage of up to 8 m²
- Wipes with tear-resistant honeycomb structure
- With fresh scent

Active ingredients: Quaternary ammonium compounds

Available as: Package with 15 wipes

Wipe size: 22 x 42 cm

Alcohol-based disinfection wipes

■ FD 322 premium wipes Quick-acting disinfection, FD 322 top wipes Quick-acting disinfection, FD 322 wipes Quick-acting disinfection

Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40.
- Very short reaction time
- Dries without residue
- Good material compatibility

Active ingredients: Alcohols

FD 322 premium wipes

Available as:

Dispenser box with 80 wipes

Wipe size: 14 x 19 cm

FD 322 wipes

Available as:

Dispenser box with 100 wipes

Wipe size: 20 x 30 cm

FD 322 top wipes

Available as: Package with 50 wipes

Wipe size: 17.5 x 25 cm

■ FD 333 wipes Quick-acting disinfection, FD 333 top wipes Quick-acting disinfection

Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses).
- Quick-acting disinfection effect
- Quick drying times
- Perfume-free

Active ingredients: Alcohols, quaternary ammonium compounds

FD 333 wipes

Available as:

Dispenser box with 110 wipes

Wipe size: 14 x 22 cm

FD 333 top wipes

Available as:

Package with 50 wipes

Wipe size: 17.5 x 25 cm

■ FD 333 forte premium wipes Quick-acting disinfection, FD 333 forte wipes Quick-acting disinfection

Fully virucidal ready-to-use disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, mycobactericidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), virucidal (enveloped and non-enveloped viruses).
- Quick-acting disinfection effect
- Extra-large wipes
- Fresh lemon scent

Active ingredients: Alcohols, quaternary ammonium compounds

FD 333 forte premium wipes

Available as:

Dispenser box with 80 wipes

Wipe size: 14 x 19 cm

FD 333 forte wipes

Available as:

Dispenser box with 100 wipes

Wipe size: 20 x 30 cm

■ FD 350 Disinfection wipes

Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40.
- Quick-acting disinfection effect
- Quick, streak-free drying
- 3 different fragrances: Classic, Flower, Lemon

Active ingredients: Alcohols

Available as: Dispenser box with 110 wipes

Wipe size: 14 x 22 cm



■ FD 366 sensitive wipes Disinfection of sensitive surfaces, FD 366 sensitive top wipes Disinfection of sensitive surfaces

Ready-to-use disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures



- Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses
- Quick-acting disinfection effect
- Extra-large wipes
- Perfume-free

Active ingredients: Alcohols, quaternary ammonium compounds

FD 366 sensitive wipes

Available as:

Dispenser box with 100 wipes

Wipe size: 20 x 30 cm

FD 366 sensitive top wipes

Available as:

Package with 50 wipes

Wipe size: 17.5 x 25 cm

Dry wipes

■ FD multi wipes

Universal dry wipes, can be soaked in disinfectant and cleaning agents



- All-purpose wipes can be used together with surface disinfectants
- Lint-free and tear-proof
- Excellent wetting; does not scratch surfaces
- Very economical

Available as: Dispenser box with 180 wipes

Wipe size: 20 x 30 cm

Our recommendation: FD 300, FD 312, FD 322, FD 333, FD 333 forte, FD 366 sensitive, wall bracket for FD multi wipes box

■ FD multi wipes compact

Universal dry wipes, can be soaked in disinfectant and cleaning agents



- All-purpose wipes can be used together with surface disinfectants
- Lint-free and tear-proof
- Excellent wetting; does not scratch surfaces
- Compact box, extra-large wipes

Available as: dispenser box with 60 wipes

Wipe size: 29 x 30 cm

Our recommendation: FD 300, FD 312, FD 322, FD 333, FD 333 forte, FD 366 sensitive

■ FD polishing wipe

Microfibre cloth with special fibres for gentle cleaning and polishing of sensitive surfaces



- Particularly thorough thanks to the fine fibre material
- Fast, streak-free cleaning
- Antistatic properties remove small particles of dirt and dust
- Suitable for washing at temperatures of up to 60°C

Available as: Individual wipe

Wipe size: 32 x 32 cm

Skin and hands



Disinfection and cleaning Hand disinfection

■ HD 410* Hand disinfection

Alcoholic rub-in preparation
for hygienic and surgical hand
disinfection



- **Hygienic hand disinfection:** bactericidal (incl. tuberculocidal / mycobactericidal), yeasticidal and effective against enveloped viruses (limited virucidal) as well as against rotaviruses, adenoviruses and noroviruses (limited virucidal PLUS).
- **Surgical hand disinfection:** bactericidal and yeasticidal.
- Extremely fast: hygienic hand disinfection in 30 seconds / surgical hand disinfection in 90 seconds
- Contains skin-care and lipid-replenishing substances, very good skin compatibility – dermatologically tested.
- Perfume-free

Active ingredients: Alcohols

Available as: 500 ml, 500 ml for Hygocare Plus, 1 l, 2.5 l, 10 l

Our recommendation: Universal dispenser S 400, L 400, L 1000,

touchless T/XT 400, T/XT 1000, RX 5 M, RX 10 M, RX 5 T, dosing pumps

Hand cleaning

■ HD 435 Washing lotion

Mild washing solution for hands



- Protects and supports the natural protective barrier of the skin – with a pH value suitable for the skin
- Dermatologically tested

Available as: 500 ml, 1 l, 2.5 l

Our recommendation: Universal dispensers S 400, L 400, L 1000,

touchless T/XT 400, T/XT 1000, RX 5 M, RX 10 M, RX 5 T, dosing pumps



Dispenser system

■ Touchless XT universal dispenser

Battery-operated sensor dispenser for contactless, hygienic application of hand preparations



- Perfect hand hygiene thanks to touchless operation
- Battery operation with an acoustic warning signal that lets you know when the battery needs to be replaced
- High standard of hygiene, thanks to an anodised aluminium housing
- Autoclavable stainless steel pump or disposable pump

Dimensions (WxHxD): XT 1000/XT 400: 9.2 x 33.0 x 21.0 cm

Available as: XT 400 for 500 ml bottles, XT 1000 for 1 l bottles

Our recommendation: HD 410, HD 435 in a 500 ml or 1 l bottle, dish holder, table stand, dispenser column, disposable pumps

■ Universal dispenser RX

Manual or sensor-activated universal dispenser for application of hand products



- Manual (RX 5 M, RX 10 M) or touchless (RX 5 T)
- Housing and collector tray made of robust white plastic
- Durable plastic pump

Dimensions (WxHxD): RX 5 M: 10.2 x 27.0 x 14.2 cm,

RX 10 M: 33.4 x 18.4 x 21.8 cm, RX 5 T: 27 x 15.7 x 22.6 cm

Available as: RX 5 M and RX 5 T for 500 ml bottle or RX 10 M for 1 l bottle

Our recommendation: HD 410, HD 435 in a 500 ml bottle or 1 l bottle

■ Universal dispensers

Universal wall-mounted dispenser for application of hand products



- Extremely user-friendly
- Available with short (S) or long (L) arm lever
- Durable aluminium and stainless steel version
- Autoclavable pump or disposable pump

Dimensions (WxHxD): S 400/L 400: 8.2 x 24.5 x 8.2 cm, arm lever length:

S: 16 cm, L: 21.5 cm/L 1000: 9.2 x 29 x 9 cm, arm lever length: 22.5 cm

Available as: S 400 and L 400 for 500 ml bottles, L 1000 for 1 l bottles

Our recommendation: HD 410, HD 435 in a 500 ml or 1 l bottle, disposable pumps

■ Dispensing station for Touchless universal dispenser

Dispenser column with stable base and frame for mounting posters relating to hand disinfection



- Freestanding
- Ideal for areas in which mounting hand sanitiser dispensers to the wall is not possible

Height: 180 cm

Oval aluminium profile, silver anodised: 7 x 3 cm

Steel base plate RAL9006 powder-coated: 42.5 x 50 cm

Aluminium folding frame in portrait format DIN A3 (29.7 x 42 cm)

Special areas



Disinfection and cleaning Suction system and spittoon

■ Orotol® plus disinfection of suction systems

Foam-free liquid concentrate for simultaneous cleaning, disinfection, deodorisation and care of all suction systems and amalgam separators



- Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), as well as non-enveloped adenoviruses and noroviruses
- Removes the most persistent deposits, even in places that are difficult to access
- The solution for contaminated suction systems – provides reliable protection against infection
- Certified material compatibility – recommended by leading equipment manufacturers

Active ingredients: Quaternary ammonium compounds

Available as: 1 l, 2.5 l

Our recommendation: OroCup, MD 555 cleaner, MD 550

■ Orotol® ultra disinfection of suction systems

Powder concentrate for simultaneous cleaning, disinfection, deodorisation and care of all suction systems and amalgam separators



- Bactericidal, tuberculocidal, mycobactericidal, fungicidal, sporicidal (Clostridium difficile), virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses, as well as non-enveloped viruses such as adenoviruses, noroviruses, polioviruses, polyomaviruses SV40).
- Removes the most persistent deposits, even in places that are difficult to access
- Low application concentration (1%)
- Certified material compatibility – recommended by leading equipment manufacturers

Active ingredients: Sodium percarbonate peroxohydrate

Available as: 500 g bag

Our recommendation: OroCup, MD 555 cleaner, MD 550

■ MD 555 cleaner:

Special detergent for suction units

Foam-free cleaning concentrate for all types of suction system and amalgam separator



- Exceptionally intensive cleaning effect against deposits
- Particularly suitable for removing residues from prophylaxis powders and prophylaxis pearl products
- Certified material compatibility – recommended by leading equipment manufacturers
- Ideal in combination with Orotol® plus

Ingredients: Organic and inorganic acids, foam-free surfactants

Available as: 2.5 l

Our recommendation: OroCup, Orotol® plus, Orotol® ultra, MD 550

■ MD 550 Spittoon bowl cleaner

Foam-free, ready-to-use solution for intensive, residue-free cleaning of spittoons



- Problem-free removal of residues, such as lime scale, blood or plaque colouring agents
- Simultaneous cleaning, deodorisation and care of spittoons
- No abrasive components – the spittoon suffers no scratching
- Anti-microbial – effective against bacteria, Tb bacteria and yeast fungus

Active ingredients: Alcohols

Available as: 750 ml

Our recommendation: Orotol® plus, Orotol® ultra, MD 555 cleaner

■ OroCup care system

User-friendly, closed care system for disinfection and cleaning of suction systems and amalgam separators



- Creates an air-water mix for optimal wetting and effective cleaning of the suction system
- Efficient work – one system for application and suctioning of the ready-to-use solutions
- Connecting adapter for all suction hoses

Available as: 2 litre care system

Our recommendation: Orotol® plus, Orotol® ultra, MD 555 cleaner



Impressions and dentistry work

■ MD 530 Cement remover and prosthesis cleaner

Ready-to-use solution for removal of cement, calculus, smoking stains, oxide layers etc. from dentures, bridges, crowns, spatulas, glass etc.



- Intensive cleaning effect and residue-free removal of acid-soluble plaque
- For direct use, outstanding material compatibility
- Particularly fast and effective, thanks to a special formulation
- Contains corrosion inhibitors
- Also suitable for use in an ultrasonic bath

Ingredients: Inorganic acids, surfactants

Available as: 2.5 l

Our recommendation: Hygobox

■ MD 535 plaster remover

Ready-to-use solution for removal of plaster and alginates from impression trays and other aids



- Neutral special cleaner with particularly fast action
- Outstanding cleaning power for residue-free removal of plaster and alginate
- With activity indicator
- Also suitable for use in an ultrasonic bath

Ingredients: Alkaline cleaning components and complexing agents

Available as: 2.5 l

Our recommendation: Hygobox

■ MD 520 Impression disinfection

Fully virucidal ready-to-use solution for effective and gentle disinfection and cleaning of impressions, prosthetic work, etc.



- Bactericidal, tuberculocidal, yeasticidal, virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses, as well as non-enveloped viruses such as adenoviruses, polyomaviruses SV40, noroviruses and polioviruses)
- For use with alginates, silicones, polyether rubber, hydrocolloids, polysulphides etc.
- Increased accuracy in the plaster master model, thanks to effective cleaning effect
- No adverse effects on dimensional stability or plaster compatibility

Active ingredients: Quaternary ammonium compounds, glutaraldehyde

Available as: 2.5 l, 10 l

Our recommendation: Disinfection box, Hygobox, Hygojet

■ MD 520 Disinfection box

Container for immersion disinfection and cleaning with MD 520.



- Contains a special suspension device in the lid for max. 2 impression trays
- With practical occupied display

Dimensions: Height 19 cm, Ø: 11.7 cm

Capacity: Maximum 1 litre

Unit

■ Hygojet

Disinfection and cleaning sheath for impressions and prosthetic work



- The closed system prevents cross-infections
- Economical to use – with just 10 ml MD 520 per impression
- Extremely effective – for 6 impressions in 12 minutes
- Tested disinfection process – proven by experts

Dimensions (WxHxD): 56 x 34 x 32 cm

Our recommendation: MD 520

Spectrums of activity and applications

	Disinfection	Cleaning	VAH listed	IHO list of disinfectants	Tested in accordance with EU standards	Bactericidal	Tuberculocidal	Yeasticidal	Fungicidal	Limited virucidal	Limited spectrum of virucidal activity	Virucidal	Concentration / application time in minutes	No. of days product works for
Instruments														
ID 212 Instrument disinfection	■	■	■	■	■	■	■	■	■	■			2%/5 ⁴	7
ID 212 forte Instrument disinfection	■	■	■	■	■	■	■	■	■	■			2%/5 ⁴	7
ID 213 Instrument disinfection	■	■	■	■	■	■	■	■	■	■			2%/5 ⁴	14
ID 215 Enzymatic instrument cleaner		■											1%/2%/30/15	
ID 220 Bur disinfection	■	■	■	■	■	■	■	■	■	■		■	100%/1 ⁴	7
Surfaces														
FD 300 Surface disinfection	■	■	■	■	■	■	■	■	■	■	■	■	1%/2%/15/5 ⁴	
FD 312 Surface disinfection ¹	■	■	■	■	■	■		■		■	■		1%/15 ⁴	
FD 322 Quick-acting disinfection	■		■	■	■	■	■	■	■	■	■		100%/1 ⁴	
FD 333 Quick-acting disinfection	■		■	■	■	■	■	■	■	■	■		100%/1 ⁴	
FD 333 forte Quick-acting disinfection	■		■	■	■	■	■	■	■	■	■	■	100%/1 ⁴	
FD 366 sensitive Disinfection of sensitive surfaces	■		■	■	■	■	■	■		■			100%/1 ⁴	
FD 312 wet wipes Surface disinfection ¹	■	■	■	■	■	■		■		■			Ready-to-use/15 ⁴	
FD 322 premium wipes for quick-acting disinfection FD 322 top wipes Quick-acting disinfection FD 322 wipes for quick-acting disinfection	■		■	■	■	■	■	■	■	■	■		Ready-to-use/1 ⁴	
FD 333 wipes for quick-acting disinfection FD 333 top wipes Quick-acting disinfection	■		■	■	■	■	■	■	■	■	■		Ready-to-use/1 ⁴	
FD 333 forte premium wipes Quick-acting disinfection FD 333 forte wipes for quick-acting disinfection	■		■	■	■	■	■	■	■	■	■	■	Ready-to-use/1 ⁴	
FD 360 Cleaning and care of vinyl upholstery		■											Ready-to-use	
FD 350 Disinfection wipes	■		■	■	■	■	■	■	■	■	■		Ready-to-use/1 ⁴	
FD 366 sensitive wipes FD 366 sensitive top wipes for disinfection of sensitive surfaces	■		■	■	■	■	■	■		■			Ready-to-use/1 ⁴	
Skin and hands														
HD 410 Hand disinfectant ^{1,2}	■		■	■	■	■	■	■	■	■	■		100%/0.5/ 1.5 ^{2,4}	
HD 435 Washing lotion		■											100%	
Special areas														
Orotol® plus disinfection of suction systems	■	■	■	■	■	■		■		■			2%/5 ⁴	
Orotol® ultra disinfection of suction systems	■	■	■	■	■	■	■	■	■	■		■	1%/5 ⁴	
MD 555 cleaner Special detergent for suction systems		■											5%/30-120	
MD 550 Spittoon bowl cleaner		■		■	■	■		■					100%/5 ⁴	
MD 520 Impression disinfection	■	■	■	■	■	■	■	■		■		■	100%/5 ^{3,4}	7
MD 530 Cement remover / cleaner of prostheses		■											100% ³	
MD 535 plaster remover		■											100%	

¹ Use disinfectants safely. Always read the label and product information before use. ² Hygienic hand disinfection/surgical hand disinfection. ³ In the Hygogel. ⁴ Bactericidal/yeasticidal. More detailed information can be found in our product information available at: www.duerrdental.com