

# Product information

## FD 312 wet wipes surface disinfection



### FD 312 wet wipes at a glance

- Ready-to-use wet wipes for cleaning and disinfection of floors, furnishings, etc. in accordance with the EU Biocidal Products Regulation (EU BPR).
- Broad spectrum of activity: bactericidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses).
- Tested in 4-field test in accordance with VAH and EN 16615.
- FD 312 is included in the VAH list.
- Good material compatibility and pleasant fragrance.
- Extra large wipes with non-tear mesh structure 220 x 420 mm.
- FD 312 wet wipes are practical and easy to use.
- Aldehyde-free – the mechanism of action is based on quaternary ammonium cations.

### Properties

FD 312 wet wipes from Dürr System Hygiene are ready-to-use wipes for cleaning and disinfection of floors, furnishings, etc. in accordance with the EU Biocidal Products Regulation (EU BPR).

### Product composition

FD 312 wet wipes contain FD 312. 100 g active substance solution contain 65 mg alkyl-benzyl-dimethylammonium chloride, nonionic surfactants, alkaline cleaning components, complexing agents, and excipients in aqueous solution.

### Antimicrobial effect

Bactericidal, yeasticidal, virucidal activity against enveloped viruses (enveloped viruses incl. HBV, HCV, HIV and coronaviruses). FD 312 is VAH-listed. Tested under high load in accordance with EN 13727, EN 13624, EN 16615.

### Directions for use

**PT 2.** Remove FD 312 wet wipes individually. Wipe the surfaces and objects to be disinfected and allow the product to take effect. Disinfect visibly contaminated surfaces with FD 312 wet wipes and remove contamination. Disinfect cleaned surfaces again. Make sure that wetting is complete throughout the entire action time. Close the flow pack immediately after removing a wipe. In accordance with VAH, the exposure time is 15 minutes (bactericidal / yeasticidal), and inactivation of enveloped viruses takes place within 60 minutes. One wipe covers approx. 6 – 8 m<sup>2</sup>.

### Environmental impact

The packaging consists of a composite material and can be thermally recycled. Dispose of used wipes with the domestic waste. Do not flush. Dispose of larger quantities in conformity with official requirements, dispose of at the hazardous or special waste collection point, as necessary. Recycle packaging. For further information see safety data sheet.

### Shelf life

Do not use the product after the expiry date printed on the packaging.

### Package size

15 wipes in the pouch

### Storage

Store in a cool, dry place in the original container. Store in a well-ventilated place. Keep container tightly closed.

### General notes

Always follow the manufacturer's instructions in the case of sur-

# Product information

## FD 312 wet wipes surface disinfection

faces made from unknown materials and/or test for compatibility in an inconspicuous place. Changes in product colour/odour may occur, especially when stored in sunlight. However, these changes in colour and/or odour have no impact on the disinfection efficacy of the product. Use exclusively by qualified personnel.

### Hazard warnings

FD 312 wet wipes are classified in accordance with the CLP Regulation and are not required to be labelled. For further information see safety data sheet.

pH: 10.3 - 11.3

### Independent expert opinions – in-house investigations

The expert opinions are available upon request.

**Use disinfectants carefully. Always read the label and summary of product characteristics before use.**

**(Biocide reg. no.: FD 312 wet wipes: N-21523)**

Application		Contact time
Recommendation for use for surface disinfection in accordance with VAH	Bactericidal / yeasticidal	15 min.

  

Type of activity	Tested according to	Contact time
<b>Bacteria and yeasts</b>		
Bactericidal / yeasticidal	EN 13727, EN 13624, EN 16615	15 min.
<b>Viruses</b>		
Virucidal activity against enveloped viruses (enveloped viruses including HBV, HCV, HIV as well as coronaviruses)	DVV / RKI guideline	60 min.

### Sales

Dürr Dental SE  
Höpfigheimer Str. 17  
74321 Bietigheim-Bissingen  
Germany  
www.duerredental.com  
info@duerrdental.com

### Manufacturer

orochemie GmbH + Co. KG  
Max-Planck-Straße 27  
70806 Kornwestheim  
Germany  
www.orochemie.de  
info@orochemie.com