

ID 212

Instrument disinfection



Aldehyde-free concentrate for the cleaning and disinfection of general/surgical instruments and rotating instruments that are sensitive to alkaline chemicals and alcohol in practices, laboratories and clinics in the dental sector



recommended by leading instrument manufacturers

fast-acting


economical

preserves the value of instruments

certified material compatibility


Exposure time

Recommendation

 immersion bath
5 min.
5 min.
ultrasonic

ideal for use in ultrasonic devices

Spectrum of activity

 **High level of protection against infection:** Bactericidal, yeasticidal, tuberculocidal in an ultrasonic device, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses, as well as non-enveloped adenoviruses and noroviruses)



Key features

