

ID 213

Instrument disinfection



Aldehyde-free concentrate for the cleaning and disinfection of general/surgical instruments and rotating instruments that are sensitive to alkaline chemicals and alcohol in practices, laboratories and clinics in the dental sector



wide spectrum of action

fast-acting

economical
(2.5 l concentrate = 125 l ready-to-use solution)

high corrosion protection

certified material compatibility – recommended by leading instrument manufacturers

Exposure time

Recommendation



immersion bath

5 min.

2 min.

ultrasonic

ideal for use in ultrasonic devices

Spectrum of activity



High level of protection against infection:

Bactericidal, tuberculocidal, mycobactericidal in an ultrasonic device, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses and polyomaviruses SV40.



Key features

14 days

service life with low application concentration

a single concentrate for all disinfection and cleaning processes

optimum material compatibility thanks to special corrosion inhibitors

