

ID 220

Bur disinfection



Aldehyde-free, ready-to-use solution for cleaning and disinfection of rotary instruments in practices, laboratories and clinics in the dental sector



ready-to-use

particularly fast-acting

exceptional cleaning performance

high corrosion protection

certified material compatibility

Exposure time



immersion bath

1 min.

30 min.

ultrasonic

Recommendation

dissolves **even the most stubborn contamination** on drills and cutters made from alkaline-resistant steel and stainless steel, carbide.

Spectrum of activity



High level of protection against infection: Bactericidal, mycobactericidal, tuberculocidal, yeasticidal, fungicidal, virucidal (enveloped viruses incl. HBV, HCV, HIV and corona-viruses, as well as non-enveloped viruses such as adeno-viruses, noroviruses, polioviruses, polyomaviruses SV40).



Certified in accordance with the DVV/RKI Directive as well as EN 14476 and EN 17111

Key features

100 %

virucidal

exceptional cleaning power, particularly for microstructures in drills and cutters

excellent material compatibility thanks to special corrosion inhibitors

