

MD 520

Impression disinfection



Highly effective formaldehyde-free, ready-to-use solution for cleaning and disinfecting dental impressions, impression trays, prosthetic work etc., for immersion disinfection in the Dürr Dental Hygobox/disinfection box.



fully virucidal

for universal use


alcohol-free

improved precision of the master cast thanks to excellent cleaning performance

preservation of dimensional stability / plaster compatibility


Exposure time

Recommendation

 just **5 min.**
VAH-compliant

can be used for impression materials, prosthetic work, orthodontic appliances, drilling templates for implant systems etc.

Spectrum of activity

 **High level of protection against infection:** Bactericidal, tuberculocidal, yeasticidal, virucidal (enveloped viruses incl. HBV, HCV, HIV, coronaviruses, as well as non-enveloped viruses such as adenoviruses, polyomaviruses SV 40, polioviruses, noroviruses).



Certified in accordance with the DVV/RKI Directive as well as EN 14476 and EN 17111

Key features

