

# MD 555 cleaner

## Special detergent for suction systems

### Cleaning of suction systems



Non-foaming special cleaner for dental suction systems including drainage lines for applications where access to tap water is difficult or drainage lines have an insufficient gradient



#### Exposure time

 **5 %**  
**30–120**  
**min.**  
at least 2 x weekly

#### Recommendation

dissolves  
**encrusted** deposits  
of **prophylaxis**  
**powders** and  
**pearl products**

#### Function



#### Descaling agent/intensive cleaning agent

Excellent cleaning performance against deposits, e.g. limescale caused by very hard water, stubborn salt residues etc.

preserves the value  
of your suction system

tested material  
compatibility

recommended by  
leading manufacturers  
of suction units

intensive cleaning effect

ideal for use  
in combination with  
Orotol® plus



#### Key features

**5 %**  
application concentration

intensive  
cleaning  
performance  
against deposits  
and sediment  
build-up

technically and  
hygienically  
perfect  
operation of  
the suction system  
is guaranteed

