

Orotol® plus

Disinfection of suction systems



Non-foaming liquid concentrate for simultaneous disinfection, deodorisation, cleaning and care of dental suction units, amalgam separators and spittoons



preserves the value of your suction system

tested material compatibility

recommended by leading manufacturers of suction units

long-lasting effects

ideal for use in combination with MD 555 cleaner

Exposure time



2%
5–60 min.
1–2 x daily

Recommendation

removes **stubborn deposits**, even in places that are **hard to get to**

Spectrum of activity



High protection against infection

Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), non-enveloped adenoviruses and noroviruses



Key features

only **2%** application concentration

dissolves and disinfects biofilm

technically and hygienically perfect operation of the suction system is guaranteed

